

iRemind

Remind yourself and your friends.

USER MANUAL

Version

0.0.1

Date

22 February 2005

Author

Teo Wee Kian

XCLUSIVE SOLUTIONS

1.0 Introduction

iRemind help to remind your friends and yourself of important dates, appointments or events via the Short-Message-Service (SMS) of the GSM network. The system is very user-friendly and easy to use.

1.1 System Requirements

- XDA O2 or XDA O2 II PDA.
- GSM subscription.
- iRemind software .

2.0 Using iRemind – Entry

Operating the iRemind consists of 5 simple steps.

Step 1: Select the date.

Step 2: Select the time.

Step 3: Enter the mobile numbers (use “;” for multiple numbers.).

Step 4: Enter the message.

Step 5: Click the “Remind” button.

2.1 iRemind - Entry Module




Figure 1

2.2 iRemind – Reminder List




Figure 2

Figure 2 shows the list of reminders – pending, partial, sent or failed.

- Pending – Reminder still waiting to be sent.
- Partial – Not all mobile numbers were sent successfully.
- Sent – All mobile numbers are sent successfully.
- Fails – Reminder failed to sent to all mobile numbers.

3.0 iRemind Configuration

iRemind comes with settings that user can customize to individual preferences. In the configuration module, there is also a tab for you to register the product with the product key.

3.1 Configuration




Figure 3

Figure 3 shows the first tab of the preferences form. This form allows you to perform quite a number of settings.

- Country Code – this is required to send the SMS successfully.
- Auto-Load keypad– if checked, the application will load the keypad whenever it focus an entry box.
- Save – remember to save before exiting the module.

XCLUSIVE SOLUTIONS

3.2 Product Registration




Figure 5

Once you installed you product. You have to e-mail the authentication code to your supplier. Within 5 working days, you will receive you product serial key. You will use this serial key to register the product.

Unregistered version only allows 30 reminders in the list. Once you hit 30 reminders, you will not able to add anymore new customer/patient.

4.0 Sales And Support

For sales enquiry, please e-mail to sales@xclusive-solutions.com

For support, please e-mail to support@xclusive-solutions.com

For more information, visit <http://www.xclusive-solutions.com>

You can also purchase iRemind at the follow websites:

- <http://www.handango.com>
- <http://www.pocketgear.com>